
CAHSS101GSMP (-40°C)AED cabinet 

- AED cabinet with GSM and AED sensor.
- An SMS is sent when the door is opened.
- Insulated stainless steel cabinet
- Frost resistant down to (-40 °C)
- Additional frontal heating element (-20 °C to -40 °C).
- Fan for ventilation.
- 2 x LED lights (meaning defibrillator is visible in the dark)
- Designed to be easy to open, even under conditions of heavy snow and frost.
- Future-proof! (Software can be upgraded)

FUNCTIONAL DESCRIPTION:

Once the door is opened... The BLINKER and SIREN will start, and GSM sends an SMS stating 'DOOR OPEN' to five selected mobile telephone numbers.
Take AED, (2 BIP) and then will be send an SMS 'NO AED'.
Place AED in the closet. (2 BIP) NewSMS 'AED PRESENT'
(If the door does not close again within 3 minutes, SIREN automatically stops)

Once the door is closed... The cabinet is once again ready for operation.

TECHNICAL DESCRIPTION:

AED sensor:	Registering if AED is in the cabinet. (IRED Infrared Sensor)
Door:	The door is fitted with thermopane glass and marked with international AED labels/symbols.
Light:	Blue blinker and 2 x LED lights so that the closet is visible in the dark.
Siren:	Siren in the cabinet. (100 dBA / 1m)
Fan:	The ventilation fan starts automatically when the temperature exceeds +40 °C.
Frost-proof cabinet:	Thermostatically-controlled heating element. Frost resistant down to -40 °C.
GSM/SMS:	Control Print with GSM for SMS communication.
Voltage:	Must be connected to a 230V electrical supply. Outdoor cable D8 mm (3 x 1.5 mm ²)
Battery backup:	12V battery backup activates in the event of power failure.
Maintenance:	A new 12V battery should be installed every three years.
External dimensions:	B x H x D = 44 x 46 x 27 cm (incl. door handles)
Internal dimensions:	B x H x D = 33 x 38 x 20 cm
Weight:	Cross weight incl. packaging 18,5 Kg.

DATA SHEET:

Connection voltage 230V. (Phase, Zero and Earth) Outdoor cable D8 mm (3 x 1.5 mm²)
Power consumption in normal operation = 40 mA (10.2 Wh)
Max. power consumption incl. Heater = 500 mA (116 Wh)



MAINTENANCE / CLEANING:

We recommend cleaning the outer surface of the cabinet at least 1-2 times per year.
Wipe using detergent and a damp cloth or brush and rinse with clean water. Finally wipe down the cabinet with a dry cloth.
If the surface of the cabinet becomes tarnished, this will be due to deposits which can easily be cleaned off.
For difficult deposits, use a scouring pad. **Never use wire wool.**

INSTALLATION - SIM card:

Remove the 2 Allen screws on the cover plate with the LED light slide the sheet out.

Install 1 pc. standard SIM card.

SIM card LED light
LED light
(Connection to telecommunications network)

Switch on the cabinet's power supply (see mounting and installation instructions)

After 30 seconds, the LED light should start to flash more slowly.
(If it blinks slowly, there is no connection to the telecommunications network)

SETTING MOBILE NUMBERS TO RECEIVE SMS:

Setting up phone numbers. Send an SMS with the following content: (+44 national code GB)

T1=+44XXXXXXXXXXXXXXXXXXXXX T2=+44XXXXXXXXXXXXXXXXXXXXX
T3=+44XXXXXXXXXXXXXXXXXXXXX T4=+44XXXXXXXXXXXXXXXXXXXXX
T5=+44XXXXXXXXXXXXXXXXXXXXX

T1=, T2=, T3= Telephone numbers that will receive the message DOOR IS OPEN.
T4=, T5= Service numbers, will receive an SMS error alert message from the closet and DOOR IS OPEN.

Receive a test SMS indicating software, temperature and whether the battery is connected: **Test**
(The cabinet will send the software version and temperature to the sender of the SMS. I.e. GSM signal OK.)
(If the battery is not connected, you will receive the message 'NO BATTERY CONNECTED')

Receive a status SMS. Send an SMS with the following content: **Status**
(The cabinet will send passwords and phone numbers to the sender.)

C1=XXXX C2=XXXX C3=XXXX C4=XXXX (Use only for enclosure locker)
T1=XXXXXXXXXXXXXXXXXXXXX T2=XXXXXXXXXXXXXXXXXXXXX
T3=XXXXXXXXXXXXXXXXXXXXX T4=XXXXXXXXXXXXXXXXXXXXX
T5=XXXXXXXXXXXXXXXXXXXXX

Receive a temperature SMS. Send an SMS with the following content: **Temp**
(The cabinet will send a message to the sender indicating the temperature, e.g. "TEMP:23". 23 C° = the current temperature in the cabinet)

ACTIVITY OF AED SENSOR:

Send a text message, **AED ON** (AED sensor enabled)
Send a text message, **AED OFF** (AED sensor not active)

SMS THE CABINET CAN SEND:

DOOR IS OPEN	If the door is opened.
NO AED	Take AED, (2 BIP) and then will be send an SMS.
AED PRESENT	Placer AED i skabet. (2 BIP) og der bliver sendt en SMS.
NO MAIN POWER RUNNING ON BATTERY POWER	230V disconnected.
MAIN POWER IS OK	230V OK.
BATTERY LOW	Battery under 11 Volt.
BATTERY DEFECTIV	Battery defectiv under 10 volt.
KEYBOARD DEFECTIVE	Keyboard must be replaced.
HEATING DEFECTIVE + TEMPERATUR	Heating element defectiv.
FAN DEFECTIVE + TEMPERATUR	Fan stopped.
TEMPERATURE HIGH	Cabinet temp 50 °C.
TEMPERATURE LOW + TEMPERATUR	Cabinet temp under 3 °C.

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